

### **Polyvinyl Alcohol**

4% aqueous: Add 40 g of polyvinyl alcohol to 1 L of hot tap water. Microwave on high for about 2 minutes; stir, and heat for additional 1–2 minute increments until dissolved. Allow solution to cool before use. (preparation of “slime”)

### **Potato Dextrose Agar**

Suspend 39 g of potato dextrose agar in 1 L of DI water. Heat to a boil while stirring constantly. Boil for 1 minute. Sterilize for 15 minutes at 121 °C (15 lbs of pressure) in an autoclave or pressure cooker. Cool to 50–55 °C and pour into sterilized culture dishes. If using for plate counts of yeasts and molds, adjust the pH to 3.5 with sterile 10% tartaric acid. (culture medium for plate counts of yeasts and molds)

Note: DI water denotes either distilled or deionized water.

ride, and 0.2 g sodium bicarbonate in 1 L of DI water. (mounting fluid and examination of blood cells)

### **Ringer’s Solution for Mammals**

Dissolve 0.42 g of potassium chloride, 9.0 g of sodium chloride, 0.24 g calcium chloride, and 0.2 g sodium bicarbonate in 1 L of DI water. (mounting fluid and examination of blood cells)

### **Rose Bengal**

1% aqueous: Dissolve 1 g of rose bengal in 50 mL DI water, then dilute to 100 mL with distilled water. Stir and filter if necessary. (biological stain)

### **Sabouraud Dextrose Agar**

Suspend 65 g of sabouraud dextrose agar in 1 L of DI water. Heat to boiling while stirring. Boil for 1 minute. Sterilize for 15

chlorinity of 19 0/00. Not for aquaria, only for technical purposes.

### **Seawater**

Dissolve 29.42 g of sodium chloride, 0.5 g of potassium chloride, 3.22 g magnesium chloride, 0.56 g sodium bromide, 1.36 g calcium sulfate, 2.4 g magnesium sulfate, 0.11 g calcium carbonate, 0.003 g ferric oxide in 1 L DI water. Not for aquaria, only for technical purposes.

### **Schiff’s Reagent**

Dissolve 0.5 g of fuchsin in 500 mL of DI water. Decolorize solution by passing sulfur dioxide gas through the solution, or add 9 g of sodium bisulfite and 20 mL of 2 M hydrochloric acid to the fuchsin solution. (test for aldehydes)

*RECIPES continued on next page.*

## **Laboratory Solutions for the Science Classroom**

***The Teacher’s Handbook to Solution Preparation***



Recipes cover biology (culture media and biological stains), chemistry and physical science solutions used in high schools. If you happen to come across a solution for which a recipe is not included, then the explanation section of the book will guide you through the steps of determining the correct procedure for