




IgG Subclass 4

Test Code

902082  

 Test information shown is for Service Area Sonora Quest Laboratories – Tempe AZ [CHANGE](#)

Clinical Significance

IgG Subclass 4

Test Resources

None found for this test

Please visit our [CLINICAL EDUCATION CENTER](#) to stay informed on any future publications, webinars, or other education opportunities.

Test Details

Methodology

Nephelometry (NEPH)

Preferred Specimen(s)

Serum (Gel Barrier) 1 mL serum from a serum separator tube (SST) (0.3 mL min). Overnight fasting is preferred.

Minimum Volume

1.00

Transport Container

Serum Separator Tube

Transport Temperature

Refrigerated

Specimen Stability

Ambient: 72 Hours Refrigerated: 3 Weeks Frozen: 6 Months

Setup Schedule

Reference ranges are provided as general guidance only. To interpret test results use the reference range in the laboratory report.

* The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.

This material contains content from LOINC® (<http://loinc.org>). The LOINC Table, LOINC Table Core are copyright © 1995-2019, Regenstrief Institute, Inc. and the Logical Observation Identifiers Names and Codes (LOINC) Committee and is available at no cost under the license at <http://loinc.org/license>.

The tests listed by specialty and category are a select group of tests offered. For a complete list of Quest Diagnostics tests, please adjust the filter options chosen, or refer to our Directory of Services.

Quest, Quest Diagnostics, the associated logo, Nichols Institute and all associated Quest Diagnostics marks are the registered trademarks of Quest Diagnostics. All third party marks— © and ™ — are the property of their respective owners. © 2000-2019 Quest Diagnostics Incorporated. All rights reserved.