

Peroxide Test

Peroxide Test

1 - 32 of 54

Related Categories: Peroxide Test Strips | Photometer | Citric Acid Strip | Oil Colorimeter | Hydrogen Peroxide

Sort Popularity ▼

View List ▼

Show 32 items ▼



Peroxide Value Photometer

Milwaukee Instruments

...the most widely used **tests** for oxidative rancidity, **peroxide** value is a measure of the concentration of **peroxides** and hydroperoxides formed in the initial stages of lipid oxidation. Milliequivalents of **peroxide** per kg of fat are measured by titration with iodide ion. **Peroxide** values are not static...

Related Products: Photometer

☒ Compare this item



Low Range Hydrogen Peroxide Test Strips

LaMotte

The Insta-Test Low Range Hydrogen **Peroxide Test Strips** is a convenient way to quickly measure **peroxide** levels in water with fast dip and read results. Range and sensitivity is 0, 1, 3, 10, 30, 50 ppm. Applications: Drinking Water, Industrial Water Boiling and Cooling, Laundry...

Related Products: Peroxide Test Strips

☒ Compare this item



PEROXIDE TEST

MilliporeSigma

Related Products: Peroxide Test Strips

☐ Compare this item

Quantofix® peroxides test sticks range, mg/L 0, 0.5, 2, 5, 10, 25

MilliporeSigma

...many **test** sticks as are required, and reseal the container immediately after use. Do not touch the **test** paper zone. Dip **test** stick briefly (approx. 1 second) into the **test** solution. After 5 seconds, compare the **test** paper zone with the

color scale. In the presence of hydrogen peroxide the test...

Related Products: Peroxide Test Strips

☐ Compare this item



Reflectoquant® Peroxide Test Method, 100-1000mg/L 50 Tests

MilliporeSigma

☐ Compare this item



Hydrogen Peroxide VACUettes Kit, Range: 0-1000 & 1200-12000 ppm

CHEMetrics

... Range: 0-1000 & 1200-12,000 ppm MDL: 200 ppm Method: Ferric Thiocyanate Kit comes in a plastic case and contains everything needed to perform 30 tests (except distilled water): Refill, Low and High Range Comparators, dilutor snapper cup, micro test tube, and instructions.

Related Products: Test Hydrogen Peroxide Kit

☐ Compare this item



Whatman Specialized Test Papers

GE Healthcare

...laboratories and in industry. Lead Acetate Test Paper Used for detecting hydrogen sulfide, this rapid qualitative test paper when wetted with distilled water can detect as little as 5 ppm of H₂S in the atmosphere or in a gas stream. Hydrogen peroxide can be detected with this paper by...

Related Products: Lead Acetate Paper

☐ Compare this item



Hydrogen Peroxide VACUettes Kit, Range: 0-100 & 120-1200 ppm

CHEMetrics

... Range: 0-100 & 120-1200 ppm MDL: 20 ppm Method: Ferric Thiocyanate Kit comes in a plastic case and contains everything needed to perform 30 tests (except distilled water): Refill, Low and High Range Comparators, dilutor snapper cup, micro test tube, and instructions.

☐ Compare this item



Hydrogen Peroxide Vacu-vials Kit, Range: 0-6.00 ppm

CHEMetrics



Hydrogen **peroxide** is a strong oxidizing agent with a variety of uses. Applications include the treating of industrial effluents and domestic waste and serving as a disinfectant in aseptic packaging. For the food and beverage industry, CHEMetrics Hydrogen **Peroxide** CHEMets® and...

☐ Compare this item

Hydrogen Peroxide SAM Kit

CHEMetrics



Hydrogen **peroxide** is a strong oxidizing agent with a variety of uses. Applications include the treating of industrial effluents and domestic waste and serving as a disinfectant in aseptic packaging. The Ferric Thiocyanate Method Reference: D. F. Boltz and J. A. Howell, eds.,...

Related Products: Hydrogen Peroxide

☐ Compare this item

Spectroquant® Hydrogen Peroxide Test Method: photometric

MilliporeSigma



UN Number: UN3316

☐ Compare this item

SafTest Peroxide Test Kit

MP Biomedicals



...**Peroxide Test Kits** are designed to measure **peroxide** content in any liquid or solid sample containing oil/fat. MP Biomedicals' **Peroxide Test Kit** contains all reagents and calibrators necessary to run samples for **peroxide** value. A control is also included in each **Peroxide Test Kit**. Each kit...

Related Products: SafTest

☐ Compare this item

MQuant™ Peroxidase, Test Method: colorimetric with test strips

MilliporeSigma



Method: colorimetric with **test** strips MQuant™ Storage class: 10 - 13 Other liquids and solids Storage: Store at +2°C to +8°C. The product quality of milk depends heavily on successful pasteurization. To verify this, the lactoperoxidase activity must be measured...

Related Products: Peroxide Test Strips

☐ Compare this item



Hydrogen Peroxide EndPoint ID® Test Kits

Thomas®

Our hydrogen **peroxide test** kit is used to verify the concentration levels of hydrogen **peroxide** used in the cleaning and sanitation industry. It is used to **test** hydrogen **peroxide** levels in sanitizer applications, wash water applications and waste water applications. Concentration levels are...

☐ Compare this item

PEROXIDE TEST

MilliporeSigma



The **peroxide test** is suitable for checking the current concentration of diluted **peroxide** solutions to ensure the disinfectant effect.

Related Products: Peroxide Test Strips

☐ Compare this item

Hydrogen Peroxide Titrets Kit

CHEMetrics



Hydrogen **peroxide** is a strong oxidizing agent with a variety of uses. Applications include the treating of industrial effluents and domestic waste and serving as a disinfectant in aseptic packaging. The Ceric Sulfate Titrimetric Method Reference: Developed by CHEMetrics, Inc. ...

☐ Compare this item

BAKER TESTSTRIPS For Peroxide (H2O2)

Avantor - J.T.Baker®



BAKER TESTSTRIPS **Test** strip for semi-quantitative determination of **Peroxide** (H₂O₂). Ideal for both lab and field use. Dip the strip and read the results from the color scale on the package. Acceptable pH of **test** solution is 5 - 7. Each container holds 100 strips.

Related Products: Peroxide Test Strips

☐ Compare this item

MQuant™ Calcium, Test Method: colorimetric with test strips and reagents

MilliporeSigma



Method: colorimetric with **test** strips and reagents 10 - 25 - 50 - 100 mg/L Ca²⁺ MQuant™ UN Number: UN3316 Storage class: 5.1B Oxidizing hazardous materials Storage: Store at +2°C to +8°C. Hazard: Corrosive The calcium **test** is suitable for rapid...

Related Products: Oxidation Test Strips

☐ Compare this item

Hydrogen Peroxide VACUettes Kit, Range: 0-25 & 30-300 ppm

CHEMetrics



... Range: 0-25 & 30-300 ppm MDL: 5 ppm Method: Ferric Thiocyanate Kit comes in a plastic case and contains everything needed to perform 30 **tests** (except distilled water): Refill, Low and High Range Comparators, dilutor snapper cup, micro **test** tube, and instructions.

☐ Compare this item

PEROXIDE TEST

MilliporeSigma



The **peroxide test** is especially suitable for monitoring prescribed cleaning agent concentrations.

Related Products: Peroxide Test Strips

☐ Compare this item

Hydrogen Peroxide CHEMets Kit, Range: 0-0.5 ppm

CHEMetrics



Hydrogen **peroxide** is a strong oxidizing agent with a variety of uses. Applications include the treating of industrial effluents and domestic waste and serving as a disinfectant in aseptic packaging. For the food and beverage industry, CHEMetrics Hydrogen **Peroxide** CHEMets® and...

Related Products: Hydrogen Peroxide

☐ Compare this item

QUANTOFIX PEROXIDES TEST STICKS 1-10MG/L

MilliporeSigma



☐ Compare this item

Hydrogen Peroxide VACUettes Kit, Range: 0-50 & 60-600 ppm

CHEMetrics



... Range: 0-50 & 60-600 ppm MDL: 10 ppm Method: Ferric Thiocyanate Kit comes in a plastic case and contains everything needed to perform 30 **tests** (except distilled water): Refill, Low and High Range Comparators, dilutor snapper cup, micro **test** tube, and instructions.

☐ Compare this item

Reflectoquant® Peroxide Test Method, 0.2-20.0mg/L:

MilliporeSigma

☐ Compare this item



Hydrogen peroxide solution contains inhibitor, 30 wt. % in H₂O, meets USP...

Honeywell

Density (g/cm³): 1.11 Molecular weight: 34.01 Molecular formula: H₂O₂ UN Number: UN2014

☐ Compare this item

Hydrogen peroxide solution

MilliporeSigma



3%, for microbiology Synonym(s): Catalase **Test** Empirical Formula: H₂O₂ Formula Weight: 34.01 Linear Formula: H₂O₂
MDL Number: MFCD00011333 Beilstein Registry Number: 3587191

☐ Compare this item

Hydrogen peroxide solution

MilliporeSigma



Contains inhibitor, 30 wt. % in H₂O, meets USP **testing** specifications Empirical Formula (Hill Notation): H₂O₂
Molecular Weight: 34.01 CAS No.: 7722-84-1 Vapor pressure: 23.3 mmHg (30 °C) Concentration: 30 wt. %...

Related Products: Hydrogen Peroxide 30%

☐ Compare this item



WaterWorks™ Low Range Peroxide Check

Industrial Test Systems

...Range **Peroxide Check** is a break through in **testing** for Hydrogen **Peroxide** (H₂O₂). Utilizing ITS' patented aperture technology, this strip allows for detection levels of 0.05, 0.3, 0.5, 1.0, 2.0, and 4.0 ppm, all in an easy dip and read procedure. Simply dip the patented aperture **test** strip into...

Related Products: Peroxide Test Strips

☐ Compare this item

Hydrogen Peroxide Vacu-vials Kit, Range: 0-3.00 ppm

CHEMetrics



Hydrogen **peroxide** is a strong oxidizing agent with a variety of uses. Applications include the treating of industrial effluents and domestic waste and serving as a disinfectant in aseptic packaging. For the food and beverage industry, CHEMetrics Hydrogen **Peroxide** CHEMets® and...

☐ Compare this item

Hydrogen Peroxide CHEMets® Kits, Ferric Thiocyanate Method

CHEMetrics



Hydrogen **peroxide** is a strong oxidizing agent with a variety of uses. Applications include the treating of industrial effluents and domestic waste and serving as a disinfectant in aseptic packaging. For the food and beverage industry, CHEMetrics Hydrogen **Peroxide** CHEMets® and...

Related Products: Sterile Sample Cups

☐ Compare this item

PERACETIC ACID TEST

MilliporeSigma



The **test** is suitable for the selective determination of peracetic acid traces in aqueous solutions and for check the absence of peracetic acid after rinsing processes, also in cases in which hydrogen **peroxide** is present.

Related Products: Peracetic Acid Test Strips

☐ Compare this item

Chloride Test Kit

LaMotte



Contents: 30 mL Hydrogen **Peroxide** 3% 30 mL Chloride Reagent A 30 mL Sulfuric Acid, 0.5N 15 mL Phenolphthalein Indicator 0.5% 2 x 30 mL Silver Nitrate 0.171 N 1 Sample Tube, 5-10-25 mL w/ cap

Related Products: Chloride Test Kit

☐ Compare this item

Related Categories: Peroxide Test Strips | Photometer | Citric Acid Strip | Oil Colorimeter | Hydrogen Peroxide

1 ▼

>

Site Search by  SLISYSTEMS.