

## Hemosiderin Stain, Urine


Test Code

518  

CPT Code(s)\*

83070

 Ordering Restrictions may apply. Please provide [SERVICE AREA INFORMATION](#) to find available tests you can order.

 Not offered in Quest Diagnostics Nichols Institute (IFD) – San Juan Capistrano. Please provide [SERVICE AREA INFORMATION](#) to find available tests you can order.

### Clinical Significance

Hemosiderin Stain, Urine - Hemosiderin may be seen in anemias caused by red cell destruction.

### Test Resources

None found for this test

Please visit our [CLINICAL EDUCATION CENTER](#) to stay informed on any future publications, webinars, or other education opportunities.

### Test Details

#### Methodology

Microscopy

#### Reference Range(s)

None detected

#### LOINC® Codes, Performing Laboratory

 *Service Area must be determined*

#### Preferred Specimen(s)

30 mL preserved urine with boric acid tablet **or**  
Unpreserved urine collected in a leak-proof container

#### Minimum Volume

20 mL

#### Transport Container

Plastic screw-cap urine container

#### Transport Temperature

**Preserved:** Room temperature  
**Unpreserved:** Refrigerated (cold packs)

#### Specimen Stability

Room temperature: 7 days  
Refrigerated: 7 days  
Frozen: Unacceptable

#### Setup Schedule

 *Service Area must be determined*

Reference ranges are provided as general guidance only. To interpret test results use the reference range in the laboratory report.

\* The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.

This material contains content from LOINC® (<http://loinc.org>). The LOINC Table, LOINC Table Core are copyright © 1995-2019, Regenstrief Institute, Inc. and the Logical Observation Identifiers Names and Codes (LOINC) Committee and is available at no cost under the license at <http://loinc.org/license>.

The tests listed by specialty and category are a select group of tests offered. For a complete list of Quest Diagnostics tests, please adjust the filter options chosen, or refer to our Directory of Services.

Quest, Quest Diagnostics, the associated logo, Nichols Institute and all associated Quest Diagnostics marks are the registered trademarks of Quest Diagnostics. All third party marks— © and ™ — are the property of their respective owners. © 2000-2019 Quest Diagnostics Incorporated. All rights reserved.