



SIMPLICITY IN WATER ANALYSIS

800.356.3072

[Home](#) [Products](#) [Technology](#) [Support](#) [News & Events](#) [About](#)

Hydrogen Peroxide — Vacu-vials® Instrumental Kit

[Product Information](#)[Instructions](#)[SDS](#)[Technical Information](#)

Hydrogen Peroxide — Vacu-vials® Instrumental Kit

Catalog No.: K-5543

Range: 0-6.00 ppm

Method: Ferric Thiocyanate

30 ampoules

Our 0-6.00 ppm Hydrogen Peroxide Test Kit comes in a cardboard box and contains everything needed to perform 30 tests: thirty ampoules, 25 mL sample cup, ampoule blank, and instructions.


This kit requires the use of the V-2000 or V-3000 Photometer, a SAM Photometer, or a spectrophotometer capable of accepting a 13 mm diameter round cell. Instrument sold separately.

Hydrogen Peroxide Analysis in the Presence of Peracetic Acid

Application Note: Aseptic Packaging/Dairy

Cost per test: \$1.02**Analysis Time:** < 1 minute**Accuracy Information:** [Download PDF](#)**Price: \$30.58**[See Discount Table](#)

Qty:

[ADD TO CART](#)

The Vacu-vials test procedure

[Click to view test procedure](#)[Zoom](#)[See Page from Catalog \(PDF\)](#)

Kit components or accessories that are sold separately (3)



Hydrogen Peroxide SAM Kit
Cat. No. I-2016



Ampoule Blank Pack
Cat. No. A-0023



Sample Cup, 25 mL
Cat. No. A-0013